

TECHNICAL PROTOCOL  
FOR  
**FRT-PGK-gb2-hygro-FRT**  
FRT flanked,  
Pro- and Eukaryotic  
Hygromycin Selection Cassette  
(A010)

Gene Bridges GmbH  
Im Neuenheimer Feld 584  
69120 Heidelberg, Germany  
Tel + 49 (0)6221 13708 11  
Fax + 49 (0)6221 13708 29  
Email: [contact@genebridges.com](mailto:contact@genebridges.com)  
[www.genebridges.com](http://www.genebridges.com)

## CONTENTS

### 1 Eppendorf tubes + manual

1. FRT-PGK-gb2-hygro-FRT: PCR template (50 ng/μl, 20μl)
2. This manual

**Store tube at -20°C**

### Please read

The products listed in this manual are for research purposes only. They are not designed for diagnostic or therapeutic use in humans, animals or plants. The Red<sup>®</sup>/ET<sup>®</sup> recombination technology is the intellectual property of Gene Bridges GmbH.

### Conditions of use

3.1 Purchaser will not manufacture, copy, reproduce, transmit, distribute, sell, lease, transfer, or improve upon the MATERIALS without prior written consent from GENE BRIDGES.

3.2 All MATERIALS relating Technologies shall be purchased from GENE BRIDGES or its authorized distributors. Use of any of the stated products from a source other than GENE BRIDGES will exempt GENE BRIDGES from any and all liabilities and warranties.

3.3 All MATERIALS purchased by research organizations, universities and other non-profit organizations may not be used for any commercial purpose. These MATERIALS are to be used for research purposes only. The MATERIALS may not be used to provide a commercial or non-commercial service, of any kind.

3.4 A purchase of MATERIALS by a private consumer is neither intended nor permitted.

## Short Description:

“FRT-PGK-gb2-hygro-FRT” template is designed to allow hygromycin selection in prokaryotic and eukaryotic cells.

The FRT-PGK-gb2-hygro-FRT template encodes the hygromycin resistance gene which combines a prokaryotic promoter (gb2) for expression in *E.coli* with a eukaryotic promoter (PGK) for expression in mammalian cells.

The prokaryotic promoter gb2 is a slightly modified version of the Em7 promoter; it mediates higher transcription efficiency than the normally used Tn5 promoter. The promoter of the mouse Phosphoglucokinase gene (PGK) is used as eukaryotic promoter. A synthetic polyadenylation signal terminates the hygromycin expression. The cassette is flanked by FRT sites for later excision by FLP-recombinase.


Using the provided PCR template one can easily create a FRT-PGK-gb2-hygro-FRT cassette flanked by any restriction sites to clone the cassette into the vector of choice. The restriction sites can be introduced by adding the corresponding sequence in the PCR primer. The template can easily be used to generate targeting constructs mediated by a single Red/ET Recombination step.

The “FRT-PGK-gb2-hygro-FRT template” is not linear but plasmid based (3643bp in size). Due to its R6K origin it can not replicate in most of the *E. coli* strains. The PCR product can therefore be used directly for downstream applications without any further purification.

We recommend growing the recombined cells between 4 and 24h under non-selective conditions before streaking the culture on plates conditioned with 100 µg/ml hygromycin plus the antibiotics for the selective marker of the target vector (e.g. 15 µg/ml chloramphenicol for BAC clones).

At least 20 PCR reactions can be performed using 1µl per reaction as template.

# Map: FRT-PGK-gb2-hygro-FRT cassette



NotI

1 AATTAACCTCACTAAGGGCGGCCGGAAGTTCCTATTCTCTAGAAAAGTATAGGAACCTTCATTTCTACCCGGTAGGGGAGG  
 82 CGCTTTTCCCAAGGCAGTCTGGAGCATGCGCTTTAGCAGCCCGCTGGGCACCTTGGCGCTACACAAGTGGCCCTCGGCCTCGCACACA  
 170 TTCCACATCCACCGGTAGGCGCCACCGGCTCGTTCCTTTGGTGGCCCGTCGCGCCACCTTCTACTCTCCCTAGTCAGGAAGTTTC  
 258 CCCCCCGCCCGCAGCTCGGCTCGTGCAGGACGTGACAAATGGAAGTAGCACGTTCTCACTAGTCTCGTGCAGATGGACAGCACCGCTG  
 346 AGCAATGGAAGCGGGTAGGCTTTGGGGCAGCGCCAAATAGCAGCTTTGCTCCTTCGCTTTCTGGGCTCAGAGGCTGGGAAGGGGTGG  
 434 GTCCGGGGCGGGCTCAGGGCGGGCTCAGGGCGGGCGGGCCCGAAGTCTCCGGAGGCCCGGCATTTGACAGCTTCAAAG  
 522 CGCAGTCTGCCCGCTGTTCTCCTCTCTCTCCTCCTCCGGCCCTTTCGACCTGCAGCAGCAGTGTGACAATTAATCATCGGCATAG  
 610 TATATCGGCATAGTATAATACGACAAGGTGAGGAACATAACC ATG AAA AAG CCT GAA CTC ACC GCG ACG TCT GTC  
 1 Met Lys Lys Pro Glu Leu Thr Ala Thr Ser Val  
 685 GAG AAG TTT CTG ATC GAA AAG TTC GAC AGC GTC TCC GAC CTG ATG CAG CTC TCG GAG GGC GAA GAA  
 12 Glu Lys Phe Leu Ile Glu Lys Phe Asp Ser Val Ser Asp Leu Met Gln Leu Ser Glu Gly Glu Glu  
 751 TCT CGT GCT TTC AGC TTC GAT GTA GGA GGG CGT GGA TAC GTC CTG CCG GTA AAT AGC TGC GCC GAT  
 34 Ser Arg Ala Phe Ser Phe Asp Val Gly Gly Arg Gly Tyr Val Leu Arg Val Asn Ser Cys Ala Asp  
 817 GGT TTC TAC AAA GAT CGT TAT GTT TAT CCG CAC TTT GCA TCG GCC GCG CTC CCG ATT CCG GAA GTG  
 56 Gly Phe Tyr Lys Asp Arg Tyr Val Tyr Arg His Phe Ala Ser Ala Ala Leu Pro Ile Pro Glu Val  
 883 CTT GAC ATT GGG GAA TTC AGC GAG AGC CTG ACC TAT TGC ATC TCC CGC CGT GCA CAG GGT GTC ACG  
 78 Leu Asp Ile Gly Glu Phe Ser Glu Ser Leu Thr Tyr Cys Ile Ser Arg Arg Ala Gln Gly Val Thr  
 949 TTG CAA GAC CTG CCT GAA ACC GAA CTG CCC GCT GTT CTG CAG CCG GTC GCG GAG GCC ATG GAT GCG  
 100 Leu Gln Asp Leu Pro Glu Thr Glu Leu Pro Ala Val Leu Gln Pro Val Ala Glu Ala Met Asp Ala  
 1015 ATC GCT GCG GCC GAT CTT AGC CAG AGC AGC GGG TTC GGC CCA TTC GGA CCG CAA GGA ATC GGT CAA  
 122 Ile Ala Ala Ala Asp Leu Ser Gln Thr Ser Gly Phe Gly Pro Phe Gly Pro Gln Gly Ile Gly Gln  
 1081 TAC ACT ACA TGG CGT GAT TTC ATA TGC GCG AIT GCT GAT CCC CAT GTG TAT CAC TGG CAA ACT GTG  
 144 Tyr Thr Thr Trp Arg Asp Phe Ile Cys Ala Ile Ala Asp Pro His Val Tyr His Trp Gln Thr Val  
 1147 ATG GAC GAC ACC GTC AGT GCG TCC GTC GCG CAG GCT CTC GAT GAG CTG ATG CTT TGG GCC GAG GAC  
 166 Met Asp Asp Thr Val Ser Ala Ser Val Ala Gln Ala Leu Asp Glu Leu Met Leu Trp Ala Glu Asp  
 1213 TGC CCC GAA GTC CCG CAC CTC GTG CAC GCG GAT TTC GGC TCC AAC AAT GTC CTG ACG GAC AAT GGC  
 188 Cys Pro Glu Val Arg His Leu Val His Ala Asp Phe Gly Ser Asn Asn Val Leu Thr Asp Asn Gly  
 1279 CGC ATA ACA GCG GTC ATT GAC TGG AGC GAG GCG ATG TTC GGG GAT TCC CAA TAC GAG GTC GCC AAC  
 210 Arg Ile Thr Ala Val Ile Asp Trp Ser Glu Ala Met Phe Gly Asp Ser Gln Tyr Glu Val Ala Asn  
 1345 ATC TTC TTC TGG AGG CCG TGG TTG GCT TGT ATG GAG CAG CAG ACG CCG TAC TTC GAG CCG AGG CAT  
 232 Ile Phe Phe Trp Arg Pro Trp Leu Ala Cys Met Glu Gln Gln Thr Arg Tyr Phe Glu Arg Arg His  
 1411 CCG GAG CTT GCA GGA TCG CCG CCG CTC CCG GCG TAT ATG CTC CCG AIT GGT CTT GAC CAA CTC TAT  
 254 Pro Glu Leu Ala Gly Ser Pro Arg Leu Arg Ala Tyr Met Leu Arg Ile Gly Leu Asp Gln Leu Tyr  
 1477 CAG AGC TTG GTT GAC GGC AAT TTC GAT GAT GCA GCT TGG GCG CAG GGT CGA TGC GAC GCA ATC GTC  
 276 Gln Ser Leu Val Asp Gly Asn Phe Asp Asp Ala Ala Trp Ala Gln Gly Arg Cys Asp Ala Ile Val  
 1543 CGA TCC GGA GCC GGG ACT GTC GGG CGT ACA CAA ATC GCC CCG AGA AGC GCG GCC GTC TGG ACC GAT  
 298 Arg Ser Gly Ala Gly Thr Val Gly Arg Thr Gln Ile Ala Arg Arg Ser Ala Ala Val Trp Thr Asp  
 1609 GGC TGT GTA GAA GTA CTT GCC GAT AGT GGA AAC CGA CCG AGC ACT CGT CCG AGG GCA AAG GAA  
 320 Gly Cys Val Glu Val Leu Ala Asp Ser Gly Asn Arg Arg Pro Ser Thr Arg Pro Arg Ala Lys Glu  
 1675 TAG GTTTCCTGCCACAGTCTGAGAGCTCCCTGGCGAATTTCGGTACCAATAAAAGAGCTTTATTTTTCATGATCTGTGTGTTGGTTTT  
 XhoI  
 1762 TGTGTGGCGCGCGAAGTTCCTATTCTCTAGAAAAGTATAGGAACCTTCCTCGAGCCCTATAGTGAGTGTATTAA